



CC14C Polyclonal Antibody

Catalog No	YP-Ab-07250
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	CDC14C CDC14B2
Immunogen	Synthesized peptide derived from human protein . at AA range: 450-530
Specificity	CC14C Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CC14C;CDC14C CDC14B2;Dual specificity protein phosphatase CDC14C (EC 3.1.3.16) (EC 3.1.3.48) (CDC14 cell division cycle 14 homolog C)
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass membrane protein . Retains its endoplasmic reticulum localization during mitosis. .
Tissue Specificity	
Function	
Background	Dual-specificity phosphatase. Preferentially dephosphorylates proteins modified by proline-directed kinases (By similarity).
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

