



PARPT Polyclonal Antibody

Catalog No	YP-Ab-07254
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	TIPARP PARP7
Immunogen	Synthesized peptide derived from human protein . at AA range: 450-530
Specificity	PARPT Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PARPT;TIPARP PARP7;TCDD-inducible poly [ADP-ribose] polymerase (EC 2.4.2.30) (ADP-ribosyltransferase diphtheria toxin-like 14) (ARTD14) (Poly [ADP-ribose] polymerase 7) (PARP-7)
Observed Band	72kD
Calculated Molecular Weight	72kD
Cell Pathway	Nucleus .
Tissue Specificity	Brain,Endometrial tumor,Testis,
Function	catalytic activity:NAD(+) + (ADP-D-ribosyl)(n)-acceptor = nicotinamide + (ADP-D-ribosyl)(n+1)-acceptor.,function:Poly [ADP-ribose] polymerase using NAD(+) as a substrate to transfer ADP-ribose onto glutamic acid residues of a protein acceptor; repeated rounds of ADP-ribosylation leads to the formation of poly(ADPribose) chains on the protein, thereby altering the function of the target protein. May play a role in the adaptative response to chemical exposure (TCDD) and thereby mediates certain effects of the chemicals.,similarity:Contains 1 C3H1-type zinc finger.,similarity:Contains 1 PARP catalytic domain.,similarity:Contains 1 WWE domain.,

Background

This gene encodes a member of the poly(ADP-ribose) polymerase superfamily. Studies of the mouse ortholog have shown that the encoded protein catalyzes histone poly(ADP-ribosyl)ation and may be involved in T-cell function. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2010],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images