



AXDN1 Polyclonal Antibody

Catalog No	YP-Ab-07259
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	AXDND1 C1orf125
Immunogen	Synthesized peptide derived from human protein . at AA range: 400-480
Specificity	AXDN1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AXDN1;AXDND1 C1orf125;Axonemal dynein light chain domain-containing protein 1
Observed Band	111kD
Calculated Molecular Weight	111kD
Cell Pathway	
Tissue Specificity	Testis,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

