



PLAP Polyclonal Antibody

Catalog No	YP-Ab-07305
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	PLAA PLAP
Immunogen	Synthesized peptide derived from human protein . at AA range: 531-580
Specificity	PLAP Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ALP;ALPP;EC:3.1.3.1;PALP;PLAP-1;PLAP;PLAA PLAP;Phospholipase A-2-activating protein (PLA2P) (PLAP)
Observed Band	66 kDa
Calculated Molecular Weight	58 kDa
Cell Pathway	Nucleus . Cytoplasm . Cell junction, synapse . Recruited to damaged lysosomes decorated with K48-linked ubiquitin chains. .
Tissue Specificity	Fetal brain,Kidney,Testis,Uterus,
Function	function:Plays an important role in the regulation of specific inflammatory disease processes.,sequence caution:Translated as Gln.,similarity:Belongs to the WD repeat PLAP family.,similarity:Contains 1 PFU domain.,similarity:Contains 1 PUL domain.,similarity:Contains 7 WD repeats.,
Background	function:Plays an important role in the regulation of specific inflammatory disease processes.,sequence caution:Translated as Gln.,similarity:Belongs to the WD repeat PLAP family.,similarity:Contains 1 PFU domain.,similarity:Contains 1 PUL domain.,similarity:Contains 7 WD repeats.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images