



F176A Polyclonal Antibody


Catalog No	YP-Ab-07331
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	FAM176A TMEM166 SP24
Immunogen	Synthesized peptide derived from human protein . at AA range: 91-140
Specificity	F176A Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F176A;FAM176A TMEM166 SP24;Protein FAM176A (Transmembrane protein 166)
Observed Band	16kD
Calculated Molecular Weight	16kD
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass membrane protein . Lysosome membrane ; Single-pass membrane protein .
Tissue Specificity	Expressed in lung, kidney, liver, pancreas, placenta, but not in heart and skeletal muscle.
Function	function:Acts as a regulator of programmed cell death, mediating both autophagy and apoptosis.,similarity:Belongs to the FAM176 family.,tissue specificity:Expressed in lung, kidney, liver, pancreas, placenta, but not in heart and skeletal muscle.,
Background	function:Acts as a regulator of programmed cell death, mediating both autophagy and apoptosis.,similarity:Belongs to the FAM176 family.,tissue specificity:Expressed in lung, kidney, liver, pancreas, placenta, but not in heart and skeletal muscle.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4°over night