




RN185 Polyclonal Antibody

Catalog No	YP-Ab-07343
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	RNF185
Immunogen	Synthesized peptide derived from human protein . at AA range: 121-170
Specificity	RN185 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RN185_HUMAN;RNF185;RING finger protein 185;2.3.2.27;RN185_MOUSE;RN185;E3 ubiquitin-protein ligase RNF185 (EC 6.3.2.-) (RING finger protein 185)
Observed Band	21kD
Calculated Molecular Weight	21kD
Cell Pathway	Mitochondrion outer membrane ; Multi-pass membrane protein . Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	Ubiquitously expressed.
Function	similarity:Contains 1 RING-type zinc finger.,
Background	similarity:Contains 1 RING-type zinc finger.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from K562 cells,
primary antibody was diluted at 1:1000, 4° over night