




DRAM1 Polyclonal Antibody

Catalog No	YP-Ab-07344
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	DRAM1 DRAM
Immunogen	Synthesized peptide derived from human protein . at AA range: 151-200
Specificity	DRAM1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DRAM1;DNA damage-regulated autophagy modulator protein 1;Damage-regulated autophagy modulator;DRAM1 DRAM;DNA damage-regulated autophagy modulator protein 1 (Damage-regulated autophagy modulator)
Observed Band	26kD
Calculated Molecular Weight	26kD
Cell Pathway	Lysosome membrane ; Multi-pass membrane protein .
Tissue Specificity	
Function	autophagy, apoptosis, cell death, programmed cell death, death,
Background	This gene is regulated as part of the p53 tumor suppressor pathway. The gene encodes a lysosomal membrane protein that is required for the induction of autophagy by the pathway. Decreased transcriptional expression of this gene is associated with various tumors. This gene has a pseudogene on chromosome 4. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4° over night