



NPB Polyclonal Antibody

Catalog No	YP-Ab-07353
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	NPB PPL7 PPNPB
Immunogen	Synthesized peptide derived from human protein . at AA range: 1-50
Specificity	NPB Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EC:1.7.-.-;Neuroglobin;Nitrite reductase;NPB;NPB PPL7 PPNPB;Neuropeptide B (Preproprotein L7) (hPPL7) [Cleaved into: Neuropeptide B-23 (NPB23) (hL7); Neuropeptide B-29 (NPB29) (hL7C)]
Observed Band	13kD
Calculated Molecular Weight	13kD
Cell Pathway	Secreted.
Tissue Specificity	Widely expressed in the central nervous system. High levels are found in substantia nigra, hypothalamus, hippocampus, spinal cord, placenta and fetal brain; lower levels are found in testis, uterus and ovary. Also detected at high levels in colorectal adenocarcinoma.
Function	function:May be involved in the regulation of feeding, neuroendocrine system, memory, learning and in the afferent pain pathway.,similarity:Belongs to the neuropeptide B/W family.,tissue specificity:Widely expressed in the central nervous system. High levels are found in substantia nigra, hypothalamus, hippocampus, spinal cord, placenta and fetal brain; lower levels are found in testis, uterus and ovary. Also detected at high levels in colorectal adenocarcinoma.,



Background

This gene encodes a member of the neuropeptide B/W family of proteins and preproprotein that is proteolytically processed to generate multiple protein products. The encoded products include neuropeptide B-23 and a C-terminally extended form, neuropeptide B-29, which are characterized by an N-terminal brominated tryptophan amino acid. Both of the encoded peptides bind with higher affinity to neuropeptide B/W (NPB/W) receptor 1 compared to the related NPB/W receptor 2. These peptides may regulate feeding, pain perception, and stress in rodents. [provided by RefSeq, Jul 2015],


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4° over night