



TRI56 Polyclonal Antibody

Catalog No	YP-Ab-07354
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	TRIM56 RNF109
Immunogen	Synthesized peptide derived from human protein . at AA range: 121-170
Specificity	TRI56 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RNF19;TRI56_HUMAN;TRIM56;RING finger protein 19;Tripartite motif-containing protein 56;2.3.2.27;TRI56;TRIM56 RNF109;E3 ubiquitin-protein ligase TRIM56 (EC 6.3.2.-) (RING finger protein 109) (Tripartite motif-containing protein 56)
Observed Band	83kD
Calculated Molecular Weight	81 kDa
Cell Pathway	Cytoplasm .
Tissue Specificity	Widely expressed (at protein level).
Function	sequence caution:Translation N-terminally extended.,sequence caution:Translation N-terminally shortened.,similarity:Belongs to the TRIM/RBCC family.,similarity:Contains 1 RING-type zinc finger.,similarity:Contains 2 B box-type zinc fingers.,
Background	sequence caution:Translation N-terminally extended.,sequence caution:Translation N-terminally shortened.,similarity:Belongs to the TRIM/RBCC family.,similarity:Contains 1 RING-type zinc finger.,similarity:Contains 2 B box-type zinc fingers.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images