



GPR97 Polyclonal Antibody

Catalog No	YP-Ab-07392
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	GPR97 PGR26
Immunogen	Synthesized peptide derived from human protein . at AA range: 210-290
Specificity	GPR97 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GPR97;PB99;PGR26;A31G24Rik;Pb1;AGRG3_HUMAN;ADGRG3;G-protein coupled receptor 97;AGRG3_MOUSE;GPR97_PGR26;Probable G-protein coupled receptor 97 (G-protein coupled receptor PGR26)
Observed Band	60kD
Calculated Molecular Weight	59 kDa
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in cultured primary dermal lymphatic endothelial cells.
Function	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 2 family. LN-TM7 subfamily.,similarity:Contains 1 GPS domain.,
Background	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 2 family. LN-TM7 subfamily.,similarity:Contains 1 GPS domain.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

