



GP182 Polyclonal Antibody


Catalog No	YP-Ab-07397
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	GPR182 ADMR
Immunogen	Synthesized peptide derived from human protein . at AA range: 190-270
Specificity	GP182 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	7TMR;ADMR;AM-R;AMR;G10D;GPR182;L1-R;gamrh;hrhAMR;G10-D;Gpcr17;Gpcr22;MB10;NOW;GP182_HUMAN;GP182_MOUSE;GP182_RAT;GP182;GPR182 ADMR;G-protein coupled receptor 182
Observed Band	44kD
Calculated Molecular Weight	44kD
Cell Pathway	Cell membrane; Multi-pass membrane protein.
Tissue Specificity	Highly expressed in heart, skeletal muscle, immune system, adrenal gland and liver.
Function	caution:Was originally (PubMed:9367907) thought to be a receptor for adrenomedullin.,function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,tissue specificity:Highly expressed in heart, skeletal muscle, immune system, adrenal gland and liver.,
Background	Adrenomedullin is a potent vasodilator peptide that exerts major effects on cardiovascular function. This gene encodes a seven-transmembrane protein that belongs to the family 1 of G-protein coupled receptors. Studies of the rat counterpart suggest that the encoded protein may function as a receptor for adrenomedullin. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night