



GP162 Polyclonal Antibody

Catalog No	YP-Ab-07404
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	GPR162 GRCA
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	GP162 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	A-2;GRCA;GP162_HUMAN;GPR162;Gene-rich cluster gene A protein; GP162_MOUSE;GP162;GPR162 GRCA;Probable G-protein coupled receptor 162 (Gene-rich cluster gene A protein)
Observed Band	63 kDa
Calculated Molecular Weight	64 kDa
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Brain,Lung,Uterus,
Function	function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 1 family.,
Background	This gene was identified upon genomic analysis of a gene-dense region at human chromosome 12p13. It appears to be mainly expressed in the brain; however, its function is not known. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images