



PACR Polyclonal Antibody

Catalog No	YP-Ab-07534
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	ADCYAP1R1
Immunogen	Synthesized peptide derived from human protein . at AA range: 80-160
Specificity	PACR Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PACR;ADCYAP1R1;Pituitary adenylate cyclase-activating polypeptide type I receptor (PACAP type I receptor) (PACAP-R-1) (PACAP-R1)
Observed Band	51kD
Calculated Molecular Weight	51kD
Cell Pathway	Cell membrane; Multi-pass membrane protein.
Tissue Specificity	Most abundant in the brain, low expression in the lung, liver, thymus, spleen, pancreas and placenta.
Function	function:This is a receptor for PACAP-27 and PACAP-38. The activity of this receptor is mediated by G proteins which activate adenylyl cyclase. May regulate the release of adrenocorticotropin, luteinizing hormone, growth hormone, prolactin, epinephrine, and catecholamine. May play a role in spermatogenesis and sperm motility. Causes smooth muscle relaxation and secretion in the gastrointestinal tract.,similarity:Belongs to the G-protein coupled receptor 2 family.,tissue specificity:Most abundant in the brain, low expression in the lung, liver, thymus, spleen, pancreas and placenta.,

Background

This gene encodes type I adenylate cyclase activating polypeptide receptor, which is a membrane-associated protein and shares significant homology with members of the glucagon/secretin receptor family. This receptor mediates diverse biological actions of adenylate cyclase activating polypeptide 1 and is positively coupled to adenylate cyclase. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified. [provided by RefSeq, Dec 2010],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images