



RGDSR Polyclonal Antibody

Catalog No	YP-Ab-07540
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	RGL4 RGR
Immunogen	Synthesized peptide derived from human protein . at AA range: 360-440
Specificity	RGDSR Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	RGDSR;RGL4 RGR;Ral-GDS-related protein (hRGR) (Ral guanine nucleotide dissociation stimulator-like 4) (RalGDS-like 4)
Observed Band	52kD
Calculated Molecular Weight	52kD
Cell Pathway	Cytoplasmic vesicle .
Tissue Specificity	PCR rescued clones,Testis,
Function	alternative products:A number of isoforms are produced,similarity:Contains 1 Ras-GEF domain.,
Background	This oncogene encodes a protein similar to guanine nucleotide exchange factor Ral guanine dissociation stimulator. Increased expression of this gene leads to translocation of the encoded protein to the cell membrane. The encoded protein can activate several pathways, including the Ras-Raf-MEK-ERK cascade. [provided by RefSeq, Jul 2016],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images