



TM7S3 Polyclonal Antibody

Catalog No	YP-Ab-07554
Isotype	IgG
Reactivity	Human;Rat;Mouse
Applications	WB;ELISA
Gene Name	TM7SF3
Immunogen	Synthesized peptide derived from human protein . at AA range: 480-560
Specificity	TM7S3 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	transmembrane 7 superfamily member 3;TM7S3;TM7SF3;Transmembrane 7 superfamily member 3 (Seven span transmembrane protein)
Observed Band	64-80 kDa
Calculated Molecular Weight	570 aa, 64 kDa
Cell Pathway	Cell membrane ; Multi-pass membrane protein .
Tissue Specificity	Widely expressed. Highly expressed in kidney and pancreas.
Function	tissue specificity:Ubiquitous. The highest expression is in kidney.,
Background	tissue specificity:Ubiquitous. The highest expression is in kidney.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

