




GPC39 Polyclonal Antibody

Catalog No	YP-Ab-07642
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	
Immunogen	Synthesized peptide derived from human protein . at AA range: 120-200
Specificity	GPC39 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GPC3;DGSX;Glypican-3 alpha subunit;Glypican-3 beta subunit;GTR2 2;GPC39; Putative G-protein coupled receptor GPCR39 (hGPCR39)
Observed Band	29kD
Calculated Molecular Weight	29kD
Cell Pathway	Cell membrane; Multi-pass membrane protein.
Tissue Specificity	
Function	
Background	GPR32P1 (G Protein-Coupled Receptor 32, Pseudogene 1) is a Pseudogene. GO annotations related to this gene include G-protein coupled receptor activity.
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from PC12 cells,
primary antibody was diluted at 1:1000, 4° over night