



CSKI2 Polyclonal Antibody

Catalog No	YP-Ab-07689
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	CASKIN2 KIAA1139
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	CSKI2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CSKI2;CASKIN2 KIAA1139;Caskin-2 (CASK-interacting protein 2)
Observed Band	132kD
Calculated Molecular Weight	132kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Brain,Epithelium,Placenta,T-cell,Testis,
Function	caution:May not bind CASK.,similarity:Contains 1 SH3 domain.,similarity:Contains 2 SAM (sterile alpha motif) domains.,similarity:Contains 6 ANK repeats.,
Background	This gene encodes a large protein that contains six ankyrin repeats, as well as a Src homology 3 (SH3) domain and two sterile alpha motif (SAM) domains, which may be involved in protein-protein interactions. The C-terminal portion of this protein is proline-rich and contains a conserved region. A related protein interacts with calcium/calmodulin-dependent serine protein kinase (CASK). Alternative splicing results in multiple transcript variants. [provided by RefSeq, Mar 2013],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images