




THYN1 Polyclonal Antibody

Catalog No	YP-Ab-07706
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	THYN1 THY28 HSPC144 MDS012 My0054
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	THYN1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HSPC144;MDS012;My0054;MY105;THY28;THYN1;THYN1 THY28 HSPC144 MDS012 My0054;Thymocyte nuclear protein 1 (Thymocyte protein Thy28)
Observed Band	24kD
Calculated Molecular Weight	24kD
Cell Pathway	Nucleus .
Tissue Specificity	Fetal brain,Spinal cord,Umbilical cord blood,Urinary bladder,
Function	PTM:Phosphorylated.,
Background	This gene encodes a protein that is highly conserved among vertebrates and plant species and may be involved in the induction of apoptosis. Alternatively spliced transcript variants encoding different isoforms have been described. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night