



TM9S1 Polyclonal Antibody

Catalog No	YP-Ab-07725
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	TM9SF1
Immunogen	Synthesized peptide derived from part region of human protein AA range: 181-230
Specificity	TM9S1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	HMP70;MP70;TM9S1_HUMAN;TM9SF1;MP70 protein family member (hMP70); transmembrane 9 superfamily member 1;TM9S1;Transmembrane 9 superfamily member 1 (MP70 protein family member) (hMP70)
Observed Band	66 kDa
Calculated Molecular Weight	66 kDa
Cell Pathway	Lysosome membrane ; Multi-pass membrane protein . Cytoplasmic vesicle, autophagosome membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in lung, pancreas, kidney, liver, placenta, skeletal muscle, heart and brain. The amount in skeletal muscle, heart and brain were considerably lower than in the other tissues.
Function	similarity:Belongs to the nonaspanin (TM9SF) family.,tissue specificity:Expressed in lung, pancreas, kidney, liver, placenta, skeletal muscle, heart and brain. The amount in skeletal muscle, heart and brain were considerably lower than in the other tissues.,
Background	similarity:Belongs to the nonaspanin (TM9SF) family.,tissue specificity:Expressed in lung, pancreas, kidney, liver, placenta, skeletal muscle, heart and brain. The amount in skeletal muscle, heart and brain were

considerably lower than in the other tissues.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images