



MVP Polyclonal Antibody

Catalog No	YP-Ab-07730
Isotype	IgG
Reactivity	Human;Mouse;Rat
Applications	WB;ELISA
Gene Name	MVP LRP
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	MVP Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	MVP;LRP;VAULT1;Major vault protein;Lung resistance-related protein;MVP LRP;Major vault protein (MVP) (Lung resistance-related protein)
Observed Band	99-104 kDa
Calculated Molecular Weight	893 aa, 99 kDa
Cell Pathway	Cytoplasm . Nucleus, nuclear pore complex . Cytoplasm, perinuclear region . 5% found in the nuclear pore complex (PubMed:15133037). Translocates from the nucleus to the cytoplasm upon EGF treatment (PubMed:16441665).
Tissue Specificity	Present in most normal tissues. Higher expression observed in epithelial cells with secretory and excretory functions, as well as in cells chronically exposed to xenobiotics, such as bronchial cells and cells lining the intestine. Overexpressed in many multidrug-resistant cancer cells.
Function	domain:MVP 3 mediates interaction with PTEN.,domain:MVP 4 mediates interaction with PARP4.,function:Required for normal vault structure. Vaults are multi-subunit structures that may act as scaffolds for proteins involved in signal transduction. Vaults may also play a role in nucleo-cytoplasmic transport. Down-regulates INFG-mediated STAT1 signaling and subsequent activation of JAK. Down-regulates SRC activity and signaling through MAP kinases.,induction:Up-regulated by IFNG.,PTM:Dephosphorylated by PTPN11.,PTM:Phosphorylated



on Tyr residues after EGF stimulation.,similarity:Contains 9 MVP (vault) repeats.,
subcellular location:5% found in the nuclear pore complex. Translocates from
the nucleus to the cytoplasm upon EGF treatment.,subunit:The vault
ribonucleoprotein particle is a huge (400 A x 670 A) cage structure of 12.9 Mda,
it consists of a dimer of half-vaults, with each half-vault co

Background

This gene encodes the major component of the vault complex. Vaults are multi-subunit ribonucleoprotein structures that may be involved in nucleo-cytoplasmic transport. The encoded protein may play a role in multiple cellular processes by regulating the MAP kinase, JAK/STAT and phosphoinositide 3-kinase/Akt signaling pathways. The encoded protein also plays a role in multidrug resistance, and expression of this gene may be a prognostic marker for several types of cancer. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, May 2012],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images