




IGKC Polyclonal Antibody

Catalog No	YP-Ab-07747
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA
Gene Name	IGKC
Immunogen	Synthesized peptide derived from part region of human protein AA range: 1-50
Specificity	IGKC Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	IgG Kappa Light Chain;IGK;IGK-;IGKC;immunoglobulin kappa locus;lg kappa chain C region
Observed Band	11kD
Calculated Molecular Weight	11kD
Cell Pathway	Secreted . Cell membrane .
Tissue Specificity	Abdominal adipose tissue,
Function	miscellaneous:The EU sequence has the INV (3) allotypic marker, Ala-45 and Val-83. The ROY sequence has the INV (1,2) allotypic marker, Ala-45 and Leu-83.,similarity:Contains 1 Ig-like (immunoglobulin-like) domain.,
Background	miscellaneous:The EU sequence has the INV (3) allotypic marker, Ala-45 and Val-83. The ROY sequence has the INV (1,2) allotypic marker, Ala-45 and Leu-83.,similarity:Contains 1 Ig-like (immunoglobulin-like) domain.,
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4°over night