



SLIT2 Polyclonal Antibody

Catalog No	YP-Ab-07762
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	SLIT2 SLIL3
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	SLIT2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Slit-2;Slit homolog 2 protein N-product;Slit homolog 2 protein C-product;Slit homolog 2 protein;Slit 2;SLIT2;SLIT2 SLIL3;Slit homolog 2 protein (Slit-2) [Cleaved into: Slit homolog 2 protein N-product;Slit homolog 2 protein C-product]
Observed Band	100 kDa, 200 kDa
Calculated Molecular Weight	170 kDa
Cell Pathway	Secreted . The C-terminal cleavage protein is more diffusible than the larger N-terminal protein that is more tightly cell associated.
Tissue Specificity	Fetal lung and kidney, and adult spinal cord. Weak expression in adult adrenal gland, thyroid, trachea and other tissues examined.
Function	domain:The leucine-rich repeat domain is sufficient for guiding both axon projection and neuronal migration, in vitro.,function:Thought to act as molecular guidance cue in cellular migration, and function appears to be mediated by interaction with roundabout homolog receptors. During neural development involved in axonal navigation at the ventral midline of the neural tube and projection of axons to different regions. SLIT1 and SLIT2 seem to be essential for midline guidance in the forebrain by acting as repulsive signal preventing inappropriate midline crossing by axons projecting from the olfactory bulb. In spinal chord development may play a role in guiding commissural axons once

they reached the floor plate by modulating the response to netrin. In vitro, silences the attractive effect of NTN1 but not its growth-stimulatory effect and silencing requires the formation of a ROBO1-DCC co

Background

This gene encodes a member of the slit family of secreted glycoproteins, which are ligands for the Robo family of immunoglobulin receptors. Slit proteins play highly conserved roles in axon guidance and neuronal migration and may also have functions during other cell migration processes including leukocyte migration. Members of the slit family are characterized by an N-terminal signal peptide, four leucine-rich repeats, nine epidermal growth factor repeats, and a C-terminal cysteine knot. Proteolytic processing of this protein gives rise to an N-terminal fragment that contains the four leucine-rich repeats and five epidermal growth factor repeats and a C-terminal fragment that contains four epidermal growth factor repeats and the cysteine knot. Both full length and cleaved proteins are secreted extracellularly and can function in axon repulsion as well as other specific processes. Alternative splici

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images