



ETS2 Polyclonal Antibody

Catalog No	YP-Ab-07776
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	ETS2
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	ETS2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ETS2;ETS2IT1;Protein C ets 2;Protein C-ets-2
Observed Band	53 kDa
Calculated Molecular Weight	469 aa, 53 kDa
Cell Pathway	Nucleus.
Tissue Specificity	Brain,Colon,Placenta,
Function	similarity:Belongs to the ETS family.,similarity:Contains 1 ETS DNA-binding domain.,similarity:Contains 1 PNT (pointed) domain.,
Background	This gene encodes a transcription factor which regulates genes involved in development and apoptosis. The encoded protein is also a protooncogene and shown to be involved in regulation of telomerase. A pseudogene of this gene is located on the X chromosome. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2012],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images