



HIC2 Polyclonal Antibody

Catalog No	YP-Ab-07808
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	HIC2 HRG22 KIAA1020 ZBTB30
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	HIC2 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Hic 2;Hic 3;HIC1-related gene on chromosome 22 protein;Hic-2;Hic-3;HIC2; HIC2 HRG22 KIAA1020 ZBTB30;Hypermethylated in cancer 2 protein (Hic-2) (HIC1-related gene on chromosome 22 protein) (Hic-3) (Zinc finger and BTB domain-containing protein 30)
Observed Band	66 kDa
Calculated Molecular Weight	615 aa, 66 kDa
Cell Pathway	Nucleus.
Tissue Specificity	Highest levels in cerebellum.
Function	function:Transcriptional repressor.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family. Hic subfamily.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 5 C2H2-type zinc fingers.,subunit:Interacts with CtBP.,tissue specificity:Highest levels in cerebellum.,
Background	function:Transcriptional repressor.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family. Hic subfamily.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 5 C2H2-type zinc fingers.,subunit:Interacts with CtBP.,tissue specificity:Highest levels in cerebellum.,

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images