



BTBD7 Polyclonal Antibody

Catalog No	YP-Ab-07857
Isotype	IgG
Reactivity	Human;Mouse
Applications	WB;ELISA
Gene Name	BTBD7 KIAA1525
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	BTBD7 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000 ELISA 1:5000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FUP1;57357E9Rik;83448M7;E13118E17Rik;mKIAA1525;BTBD7_HUMAN; BTBD7;KIAA1525;BTBD7_MOUSE;Function-unkown protein 1;BTB domain containing 7;BTB (POZ) domain containing 7;BTBD7 KIAA1525;BTB/POZ domain-containing protein 7
Observed Band	124kD
Calculated Molecular Weight	126 kDa
Cell Pathway	Nucleus .
Tissue Specificity	Amygdala,Brain,Neuroblastoma,Testis,
Function	similarity:Contains 2 BTB (POZ) domains.,
Background	similarity:Contains 2 BTB (POZ) domains.,
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images

