




BRX1 rabbit pAb

Catalog No	YP-Ab-07914
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	BRX1 BRX BXDC2
Immunogen	Synthesized peptide derived from human BRX1 AA range: 15-65
Specificity	This antibody detects endogenous levels of BRX1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.28% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BRCC1; BROVCA1; EC:2.3.2.27; IRIS; PSCP; BRX1; BRX1 BRX BXDC2; Ribosome biogenesis protein BRX1 homolog (Brix domain-containing protein 2)
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Nucleus, nucleolus.
Tissue Specificity	Cervix, Placenta,
Function	function: Required for biogenesis of the 60S ribosomal subunit., similarity: Belongs to the BXDC2/BRX1 family., similarity: Contains 1 Brix domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4° over night