




C2D4D rabbit pAb

Catalog No	YP-Ab-07915
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	C2CD4D FAM148D
Immunogen	Synthesized peptide derived from human C2D4D AA range: 254-304
Specificity	This antibody detects endogenous levels of C2D4D at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.29% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FAM148D;C2D4D_HUMAN;C2CD4D;C2 calcium dependent domain containing 4D;family with sequence similarity 148;member D;C2D4D;C2CD4D FAM148D; C2 calcium-dependent domain-containing protein 4D
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	plasma membrane,
Tissue Specificity	Brain,
Function	similarity:Contains 1 C2 domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4° over night