




CC032 rabbit pAb

Catalog No	YP-Ab-07916
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	C3orf32 FLS485
Immunogen	Synthesized peptide derived from human CC032 AA range: 273-323
Specificity	This antibody detects endogenous levels of CC032 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.30% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CC032;C3orf32 FLS485;Uncharacterized protein C3orf32
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Expressed in enterocytes of small and large intestinal mucosa (at protein level). Expressed in chromaffine and interstitial cells. Expressed in peripheral blood and gingival cells (PubMed:27680507).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from DU145 cells,
primary antibody was diluted at 1:1000, 4° over night