




TCP2L rabbit pAb

Catalog No	YP-Ab-07927
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	TCP10L2
Immunogen	Synthesized peptide derived from human TCP2L AA range: 276-326
Specificity	This antibody detects endogenous levels of TCP2L at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.41% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TCP2L;TCP10L2;Putative t-complex protein 10A homolog 2 (T-complex 10-like protein 2)
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	
Tissue Specificity	Testis,
Function	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data., similarity:Belongs to the TCP10 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 3T3 cells, primary antibody was diluted at 1:1000, 4° over night