




THUM1 rabbit pAb

Catalog No	YP-Ab-07928
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	THUMPD1
Immunogen	Synthesized peptide derived from human THUM1 AA range: 91-141
Specificity	This antibody detects endogenous levels of THUM1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.42% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	THUM1;THUMPD1;THUMP domain-containing protein 1
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	nucleoplasm,
Tissue Specificity	Endometrium,Epithelium,Liver,Lung,
Function	similarity:Belongs to the THUMPD1 family.,similarity:Contains 1 THUMP domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night