




TSK rabbit pAb

Catalog No	YP-Ab-07930
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	TSKU E2IG4 LRRC54 TSK UNQ850/PRO1788
Immunogen	Synthesized peptide derived from human TSK AA range: 288-338
Specificity	This antibody detects endogenous levels of TSK at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.44% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Kinase EMT; Interleukin-2-inducible T-cell kinase; IL-2-inducible T-cell kinase; TSK; TSKU E2IG4 LRRC54 TSK UNQ850/PRO1788; Tsukushin (Tsukushi) (E2-induced gene 4 protein) (Leucine-rich repeat-containing protein 54)
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Secreted .
Tissue Specificity	Brain, Colon, Placenta, Uterus,
Function	induction: By 17-beta-estradiol., similarity: Contains 9 LRR (leucine-rich) repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from VEC cells,
primary antibody was diluted at 1:1000, 4° over night