




ARHL1 rabbit pAb

Catalog No	YP-Ab-07931
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	ADPRHL1 ARH2
Immunogen	Synthesized peptide derived from human ARHL1 AA range: 174-224
Specificity	This antibody detects endogenous levels of ARHL1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.45% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ARHL1; ADPRHL1 ARH2; [Protein ADP-ribosylarginine] hydrolase-like protein 1 (ADP-ribosylhydrolase 2)
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Cytoplasm, myofibril, sarcomere .
Tissue Specificity	Ovary,
Function	similarity: Belongs to the ADP-ribosylglycohydrolase family.,
Background	ADP-ribosylation is a reversible posttranslational modification used to regulate protein function. ADP-ribosyltransferases (see ART1; MIM 601625) transfer ADP-ribose from NAD ⁺ to the target protein, and ADP-ribosylhydrolases, such as ADPRHL1, reverse the reaction (Glowacki et al., 2002 [PubMed 12070318]) .[supplied by OMIM, Mar 2008],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from THP-1 cells, primary antibody was diluted at 1:1000, 4° over night