



ATOH1 rabbit pAb


Catalog No	YP-Ab-07932
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ATOH1 ATH1 BHLHA14
Immunogen	Synthesized peptide derived from human ATOH1 AA range: 38-88
Specificity	This antibody detects endogenous levels of ATOH1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.46% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ATH1;bHLHa14;Class A basic helix-loop-helix protein 14;HATH1;Helix-;ATOH1; ATOH1 ATH1 BHLHA14;Protein atonal homolog 1 (Class A basic helix-loop-helix protein 14) (bHLHa14) (Helix-loop-helix protein hATH-1) (hATH1)
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Nucleus .
Tissue Specificity	PCR rescued clones,
Function	function:Activates E box-dependent transcription in collaboration with TCF3/E47, but the activity is completely antagonized by the negative regulator of neurogenesis HES1. May play a role in the differentiation of subsets of neural cells by activating E box-dependent transcription.,similarity:Contains 1 basic helix-loop-helix (bHLH) domain.,subunit:Efficient DNA binding requires dimerization with another bHLH protein.,
Background	This protein belongs to the basic helix-loop-helix (BHLH) family of transcription factors. It activates E-box dependent transcription along with E47. [provided by RefSeq, Jul 2008],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HCT116 cells, primary antibody was diluted at 1:1000, 4° over night