




CC037 rabbit pAb

Catalog No	YP-Ab-07934
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	C3orf37 DC12
Immunogen	Synthesized peptide derived from human CC037 AA range: 173-223
Specificity	This antibody detects endogenous levels of CC037 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.48% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CC037;C3orf37 DC12;UPF0361 protein C3orf37
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Chromosome . Recruited to chromatin following DNA damage (PubMed:30554877). Localizes to replication forks (PubMed:30554877). .
Tissue Specificity	Brain,Dendritic cell,Lung,Uterus,
Function	similarity:Belongs to the UPF0361 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4° over night