




F132B rabbit pAb

Catalog No	YP-Ab-07937
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM132B
Immunogen	Synthesized peptide derived from human F132B AA range: 102-152
Specificity	This antibody detects endogenous levels of F132B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.51% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F132B;FAM132B;Protein FAM132B
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Secreted .
Tissue Specificity	Cerebellum,Testis,
Function	similarity:Belongs to the FAM132 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from VEC cells,
primary antibody was diluted at 1:1000, 4° over night