




KLD8B rabbit pAb

Catalog No	YP-Ab-07942
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	KLHDC8B FP17659
Immunogen	Synthesized peptide derived from human KLD8B AA range: 120-170
Specificity	This antibody detects endogenous levels of KLD8B at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.56% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KLD8B; KLHDC8B FP17659; Kelch domain-containing protein 8B
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Cytoplasm . Midbody . In mitotic cells, concentrates in the midbody of the cytoplasmic bridge linking daughter cells as they are about to separate during cytokinesis. .
Tissue Specificity	Adrenal cortex, Brain,
Function	similarity: Contains 8 Kelch repeats.,
Background	This gene encodes a protein which forms a distinct beta-propeller protein structure of kelch domains allowing for protein-protein interactions. Mutations in this gene have been associated with Hodgkin lymphoma. [provided by RefSeq, Sep 2010],
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from PC-12 cells, primary antibody was diluted at 1:1000, 4° over night