




KT112 rabbit pAb

Catalog No	YP-Ab-07943
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	KT112 SBBI81
Immunogen	Synthesized peptide derived from human KT112 AA range: 98-148
Specificity	This antibody detects endogenous levels of KT112 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.57% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TOT4; SBBI81; KT1 12; KT112; KT112 SBBI81; Protein KT112 homolog
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	
Tissue Specificity	Skin,
Function	similarity: Belongs to the KT112 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from AD293 cells,
primary antibody was diluted at 1:1000, 4° over night