




DCD2C rabbit pAb

Catalog No	YP-Ab-07951
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	DCDC2C
Immunogen	Synthesized peptide derived from human DCD2C AA range: 266-316
Specificity	This antibody detects endogenous levels of DCD2C at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.65% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DCD2C;DCDC2C;Doublecortin domain-containing protein 2C
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Cell projection, cilium, flagellum . Cytoplasm . Detected along the length of the sperm flagellum and in the cytoplasm of the germ cells. .
Tissue Specificity	Expressed in testis and spermatozoa (at protein level).
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from PC-12 cells,
primary antibody was diluted at 1:1000, 4° over night