




IGFR1 rabbit pAb

Catalog No	YP-Ab-07955
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	IGFLR1 TMEM149 U2AF1L4
Immunogen	Synthesized peptide derived from human IGFR1 AA range: 159-209
Specificity	This antibody detects endogenous levels of IGFR1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.69% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CD64A;Fc gamma RI;Fc gamma RIA;FCG1;Fc-gamma RI;IGFR1;IGFLR1 TMEM149 U2AF1L4;IGF-like family receptor 1 (Transmembrane protein 149) (U2 small nuclear RNA auxiliary factor 1-like 4)
Observed Band	38kD
Calculated Molecular Weight	38kD
Cell Pathway	Cell membrane ; Single-pass type I membrane protein .
Tissue Specificity	Eye,Small intestine,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night