




LRIQ4 rabbit pAb

Catalog No	YP-Ab-07961
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	LRRIQ4 LRRC64
Immunogen	Synthesized peptide derived from human LRIQ4 AA range: 480-530
Specificity	This antibody detects endogenous levels of LRIQ4 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.75% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	LRIQ4;LRRIQ4 LRRC64;Leucine-rich repeat and IQ domain-containing protein 4 (Leucine-rich repeat-containing protein 64)
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	
Tissue Specificity	
Function	similarity:Contains 1 IQ domain.,similarity:Contains 20 LRR (leucine-rich) repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MDA-MB cells,
primary antibody was diluted at 1:1000, 4° over night