




ZNF34 rabbit pAb

Catalog No	YP-Ab-07968
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZNF34 KOX32
Immunogen	Synthesized peptide derived from human ZNF34 AA range: 203-253
Specificity	This antibody detects endogenous levels of ZNF34 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.82% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KOX32;zinc finger protein 34;Zinc finger protein KOX32;ZNF34;ZNF34 KOX32; Zinc finger protein 34 (Zinc finger protein KOX32)
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Nucleus .
Tissue Specificity	Blood,Brain,Lung,Substantia nigra,
Function	function:May be involved in transcriptional regulation.,sequence caution:Codon in position 56 corresponds to a stop codon in genome.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 12 C2H2-type zinc fingers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4°over night