




TM209 rabbit pAb

Catalog No	YP-Ab-07977
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	TMEM209
Immunogen	Synthesized peptide derived from human TM209 AA range: 421-471
Specificity	This antibody detects endogenous levels of TM209 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.91% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TM209; TMEM209; Transmembrane protein 209
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Brain, Teratocarcinoma, Testis,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from THP-1 cells,
primary antibody was diluted at 1:1000, 4°over night