




ZN394 rabbit pAb

Catalog No	YP-Ab-07979
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	ZNF394 ZKSCAN14
Immunogen	Synthesized peptide derived from human ZN394 AA range: 279-329
Specificity	This antibody detects endogenous levels of ZN394 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.93% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZN394;ZNF394 ZKSCAN14;Zinc finger protein 394 (Zinc finger protein with KRAB and SCAN domains 14)
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Nucleus .
Tissue Specificity	Colon,Fetal brain,Kidney,Lung,Mammary gland,Skin,
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 1 SCAN box domain.,similarity:Contains 7 C2H2-type zinc fingers.,
Background	The protein encoded by this gene is a zinc finger protein that inhibits the transcription of mitogen-activated protein kinase signaling pathways. The encoded protein may be involved in cardiac function. [provided by RefSeq, Sep 2016],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4°over night