




F222B rabbit pAb

Catalog No	YP-Ab-07981
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM222B C17orf63
Immunogen	Synthesized peptide derived from human F222B AA range: 261-311
Specificity	This antibody detects endogenous levels of F222B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.95% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F222B;FAM222B C17orf63;Protein FAM222B
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	nucleoplasm,
Tissue Specificity	Colon,Kidney epithelium,Lung,Teratocarcinoma,Trachea,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from VEC cells,
primary antibody was diluted at 1:1000, 4° over night