




CCD63 rabbit pAb

Catalog No	YP-Ab-07991
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	CCDC63
Immunogen	Synthesized peptide derived from human CCD63 AA range: 25-75
Specificity	This antibody detects endogenous levels of CCD63 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.105% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CCD63; CCDC63; Coiled-coil domain-containing protein 63
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	axoneme,
Tissue Specificity	Brain, Testis,
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night