




NXP20 rabbit pAb

Catalog No	YP-Ab-07996
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM114A1 NOXP20
Immunogen	Synthesized peptide derived from human NXP20 AA range: 135-185
Specificity	This antibody detects endogenous levels of NXP20 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.110% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	NXP20;FAM114A1 NOXP20;Protein NOXP20 (Nervous system overexpressed protein 20) (Protein FAM114A1)
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Cytoplasm .
Tissue Specificity	Colon,Kidney,Salivary gland,Trachea,
Function	function:May play a role in neuronal cell development.,similarity:Belongs to the FAM114 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SW480 cells,
primary antibody was diluted at 1:1000, 4° over night