



ZN461 rabbit pAb


Catalog No	YP-Ab-08000
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	ZNF461 GIOT1
Immunogen	Synthesized peptide derived from human ZN461 AA range: 283-333
Specificity	This antibody detects endogenous levels of ZN461 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.114% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GIOT-1;GIOT1;Gonadotropin inducible ovary transcription repressor 1; Gonadotropin inducible transcription repressor 1;Gonadotropin-inducible ovary transcription repressor 1;Zinc finger protein 461;ZN461_HUMAN;ZNF461; ZN461;ZNF461 GIOT1;Zinc finger protein 461 (Gonadotropin-inducible ovary transcription repressor 1) (GIOT-1)
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Nucleus .
Tissue Specificity	Widely expressed, with highest levels in liver, kidney, pancreas, thymus, and small intestine.
Function	function:May be involved in transcriptional regulation.,similarity:Belongs to the krueppel C2H2-type zinc-finger protein family.,similarity:Contains 1 KRAB domain.,similarity:Contains 12 C2H2-type zinc fingers.,tissue specificity:Widely expressed, with highest levels in liver, kidney, pancreas, thymus, and small intestine.,
Background	

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products ImagesWestern blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4°over night