




ZN879 rabbit pAb

Catalog No	YP-Ab-08002
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ZNF879
Immunogen	Synthesized peptide derived from human ZN879 AA range: 446-496
Specificity	This antibody detects endogenous levels of ZN879 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.116% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZN879;ZNF879;Zinc finger protein 879
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Nucleus .
Tissue Specificity	Prostate,
Function	similarity:Contains 1 KRAB domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4° over night