




ANTRL rabbit pAb

Catalog No	YP-Ab-08010
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ANTXRL
Immunogen	Synthesized peptide derived from human ANTRL AA range: 332-382
Specificity	This antibody detects endogenous levels of ANTRL at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.124% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ANTRL;ANTXRL;Anthrax toxin receptor-like
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A549 cells,
primary antibody was diluted at 1:1000, 4°over night