



SNUT2 rabbit pAb

Catalog No	YP-Ab-08014
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	USP39 CGI-21 HSPC332 PRO2855
Immunogen	Synthesized peptide derived from human SNUT2 AA range: 92-142
Specificity	This antibody detects endogenous levels of SNUT2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.128% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	65K;CGI 21;CGI-21;FLJ33136;HSPC332;Inactive ubiquitin specific peptidase 39; Inactive ubiquitin-specific peptidase 39;MGC7569;PRO2855;S. CEREVISIAE; HOMOLOG OF SAD1;SAD1;SAD1 homolog;Small nuclear ribonucleoprotein 65kDa (U4/U6.U5);snRNP ASSEMBLY DEFECTIVE 1;SnRNP assembly defective 1 homolog;SNRNP65;SNUT2_HUMAN;U4/U6.U5 tri snRNP associated 65 kDa protein;U4/U6.U5 tri snRNP associated protein 2;U4/U6.U5 tri-snRNP-associated 65 kDa protein;U4/U6.U5 tri-snRNP-associated protein 2; Ubiquitin specific peptidase 39;Ubiquitin specific protease 39;USP39;SNUT2; USP39 CGI-21 HSPC332 PRO2855;U4/U6.U5 tri-snRNP-associated protein 2 (Inactive ubiquitin-specific peptidase 39) (SAD1 homolog) (U4/U6.U5 tri-snRNP-associated 65 kDa protein) (65K)
Observed Band	60kD
Calculated Molecular Weight	60kD
Cell Pathway	Nucleus .
Tissue Specificity	Epithelium,Fetal liver,Liver,Skin,Umbilical cord blood,Uterus,
Function	function:May play a role in mRNA splicing. It is unsure if the protein really exhibits hydrolase activity. Could be a competitor of ubiquitin C-terminal hydrolases (UCHs).,similarity:Belongs to the peptidase C19 family.,

similarity:Contains 1 UBP-type zinc finger.,subunit:Part of a tri-snRNP complex.,

Background


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night